




Process to update or format current license

Prerequisites: RASWin Licensing Tool. Could be found on RTE folder where RASWin is installed a shortcut to it is also present in the RASWin application menu from Windows.



It also could be downloaded from:


<http://www.raswin.eu/raswinupdates/RASWinHasp.zip>

Proceed to next page.




Process to update your current license with a VC2.

- 1)** Execute tool RASWinHasp.exe
- 2)** In the dialog select tab **Apply License File**.
- 3)** In the bottom of the dialog press the browse file button [...] and locate the V2C file.
- 4)** Click **Apply Update**.



After that your license should be update and you can run RASWin normally. If you receive and error like this:




That means you need to install the RTE Vendor library. If this happens ask to us for the RTE installer or you can download it from our ftp at <http://www.raswin.eu/raswinupdates/VendorRTE.zip>.

When you've downloaded the zip file, you'll have to unzip it in a folder of your choice, then open and you'll see two files:

Nombre	Fecha de modifica...	Tipo	Tamaño
haspdinst.exe	08/09/2016 14:33	Aplicación	16.763 KB
update_vendor.bat	08/09/2016 14:31	Archivo por lotes ...	1 KB

Execute the **update_vendor.bat** file and wait for the process to complete, this message should appear.





Now [restart](#) the update/format license process